

CLICKIT ER

EMERGENCY EXTERNAL FIXATION SYSTEM

Universal clamp

Snap and Click Coupling

MR Conditional

Sterile



Indications

ClickIt ER is an MR Conditional emergency external fixation system for temporary stabilisation of open fractures, unstable fractures and fractures associated with soft tissue injury.

FEATURES

Universal clamp

Provisional tightening by hand

Snap and Click Coupling

Only one instrument required

Sterile

MR Conditional

ADVANTAGES

Rod/rod, rod/screw (Ø 4.0 - Ø 5.0 - Ø 6.0 mm) and rod/wire connection

Quick application

Quick application

Simple

Ready to use

Allows up to 3 Tesla magnetic resonance

Components




Universal clamp


One clamp for rod/rod, rod/screw, rod/wire connections



Rod/rod connection Ø 12 mm



Rod/screw connection Ø 12 mm (rod) and Ø 4, 5, 6 mm (screw)



Rod/wire connection Ø 12 mm (rod) and from 1.8 to 2.5 mm (wire)




Rods

Diameter 12 mm radiolucent



Multiscrew clamp

Easy positioning of 2 or 3 screws



Screws and Wires

Available in different diameter, material and thread lengths




ER hybrid connecting Post

For ring connection in hybrid configurations




"U" rods

Optimal accessory for tibial pilon and tibial plateau




Arch 5/8

For hybrid configurations



Half rings

For hybrid configurations



Radiolucent carbon ring

For hybrid configurations

Custom Kit ER




Code	Description	Q.ty
5006508ST	Universal clamp	max 6



Code	Description	Q.ty
5006504ST	Multiscrew clamp	max 2



Code	Description
5006510ST	Rod Ø 12 mm L 100 mm
5006512ST	Rod Ø 12 mm L 150 mm
5006514ST	Rod Ø 12 mm L 200 mm
5006516ST	Rod Ø 12 mm L 250 mm
5006518ST	Rod Ø 12 mm L 300 mm
5006520ST	Rod Ø 12 mm L 350 mm




Code	Description
5000592ST	Transfixation screw L 250 mm
VA4.150.40	Screw Ø 4.0 mm L 150x40 mm
VA5.150.40	Screw Ø 5.0 mm L 150x40 mm
VA6.250.40	Screw Ø 6.0 mm L 250x40 mm

Configurations


UPPER LIMB






PELVIS




LOWER LIMB










Kit ClickIt ER

ANKLE SPANNING KIT



5006600ST

Kit consisting of	Q.ty
Multiscrew clamp	1
Universal clamp	4
Rod Ø 12 mm L 200 mm	1
Rod Ø 12 mm L 350 mm	2
Transfixation screw L 250 mm	1
Screw Ø 5.0 mm L 150x40 mm	2

ANKLE SPANNING KIT PLUS



5006600STB

Kit consisting of	Q.ty
Multiscrew clamp	1
Universal clamp	6
Rod Ø 12 mm L 200 mm	2
Rod Ø 12 mm L 350 mm	2
Transfixation screw L 250 mm	1
Screw Ø 4.0 mm L 150x40 mm	1
Screw Ø 5.0 mm L 150x40 mm	2


KNEE BRIDGING KIT



5006630ST01

Kit consisting of	Q.ty
Universal clamp	6
Rod Ø 12 mm L 250 mm	1
Rod Ø 12 mm L 350 mm	2
Screw Ø 5.0 mm L 150x40 mm	2
Screw Ø 6.0 mm L 250x40 mm	2


ELBOW KIT PLUS



5006650STA

Kit consisting of	Q.ty
Universal clamp	7
Rod Ø 12 mm L 250 mm	2
Rod Ø 12 mm L 300 mm	1
Screw Ø 4.0 mm L 150x40 mm	2
Screw Ø 5.0 mm L 150x40 mm	2


ANKLE SPANNING KIT L



5006600STA

Kit consisting of	Q.ty
MultiscREW clamp	1
Universal clamp	4
Rod Ø 12 mm L 200 mm	1
Rod Ø 12 mm L 350 mm	2
Transfixation screw L 250 mm	1
Screw Ø 5.0 mm L 150x40 mm	2


TIBIA KIT



5006620ST01

Kit consisting of	Q.ty
Universal clamp	6
Rod Ø 12 mm L 150 mm	2
Rod Ø 12 mm L 200 mm	1
Screw Ø 5.0 mm L 150x40 mm	4


FEMUR KIT



5006640ST01

Kit consisting of	Q.ty
Universal clamp	6
Rod Ø 12 mm L 150 mm	1
Rod Ø 12 mm L 200 mm	2
Screw Ø 6.0 mm L 250x40 mm	4


HUMERUS KIT



5006660ST01

Kit consisting of	Q.ty
Universal clamp	6
Rod Ø 12 mm L 150 mm	3
Screw Ø 5.0 mm L 150x40 mm	4

ANKLE BRIDGING KIT DELTA PLUS



5006610ST02

Kit consisting of	Q.ty
Universal clamp	8
Rod Ø 12 mm L 200 mm	1
Rod Ø 12 mm L 250 mm	1
Rod Ø 12 mm L 300 mm	2
Transfixation screw L 250 mm	1
Screw Ø 5.0 mm L 150x40 mm	2
Screw Ø 4.0 mm L 150x40 mm	1

ClickIt ER instrument set



ClickIt ER single use instrument set




5006010ST01

Code	Description	Q.ty
5000603M	Universal T handle for screw	1
5006010M	Guide Ø 6.0 L 60 - Single use	1
5006000M	Hexagonal L-wrench 3.0 mm - Single use	1
5002074M	Hexagonal L-wrench 5.0 mm - Single use	1

5006012ST01

Code	Description	Q.ty
5000603M	Universal T handle for screw	1
5006012M	Guide Ø 6.0 L 100 - Single use	1
5006000M	Hexagonal L-wrench 3.0 mm - Single use	1
5002074M	Hexagonal L-wrench 5.0 mm - Single use	1

ClickIt ER non sterile set



Components

Code	Description	Q.ty
5006508	Universal clamp	14
5006504	ER multiscrew clamp	2
5006510	Rod Ø 12 L 100 mm	2
5006512	Rod Ø 12 L 150 mm	2
5006514	Rod Ø 12 L 200 mm	2
5006516	Rod Ø 12 L 250 mm	2
5006518	Rod Ø 12 L 300 mm	2
5006520	Rod Ø 12 L 350 mm	2
5006522	Rod Ø 12 L 400 mm	2
5.2201A.002	ClickIt ER "U" rod ER Ø 180	1
5006560M	Metal set for ClickIt ER system	1

CLICKIT ER

COMPONENTS and ACCESSORIES

Code	Description
5006508ST	Universal clamp ER
5006504ST	ER multiscrew clamp
5006506ST	Hybrid Connecting Post (ClickIt ER)
5006510ST	Rod Ø 12 L 100 mm
5006512ST	Rod Ø 12 L 150 mm
5006514ST	Rod Ø 12 L 200 mm
5006516ST	Rod Ø 12 L 250 mm
5006518ST	Rod Ø 12 L 300 mm
5006520ST	Rod Ø 12 L 350 mm
5006522ST	Rod Ø 12 L 400 mm
5006526ST	Rod Ø 12 L 500 mm
5006528ST	Rod Ø 12 L 600 mm
5.2201A.001ST	ClickIt ER “U” rod ER Ø 150
5.2201A.002ST	ClickIt ER “U” rod ER Ø 180
5.2201A.003ST	ClickIt ER “U” rod ER Ø 210
5006102ST	5/8 ring Ø 120 mm
5006104ST	5/8 ring Ø 140 mm
5006106ST	5/8 ring Ø 160 mm
5006108ST	5/8 ring Ø 180 mm
5006110ST	5/8 ring Ø 200 mm
5006112ST	5/8 ring Ø 220 mm

CLICKIT ER KIT

5006600ST	ClickIt ER Ankle Bridging kit
5006600STB	ClickIt ER Ankle Bridging kit Plus
5006610ST02	Ankle Bridging kit Delta Plus
5006620ST01	ClickIt ER Tibia Kit
5006630ST01	ClickIt ER Knee Bridging Kit
5006640ST01	ClickIt ER Femur Kit
5006650STA	ClickIt ER Elbow Kit Plus
5006660ST01	ClickIt ER Humerus Kit

INSTRUMENTS *

Code	Description	Q.ty
5006562	Complete ClickIt ER instrument tray	
Comprising:		
5000124	3.0 mm hex screwdriver	1
5000600	T-wrench for screws Ø 4.0 mm	1
5000601	T-wrench for screws Ø 5.0 mm	1
5000602	T-wrench for screws Ø 6.0 mm	1
5000604	T-handle with AO attachment	1
5000605	Hexagonal rod Ø 5.0 mm AO connection	1
5006000	3.0 mm L-shaped Allen key	1
5002074	5.0 mm L-shaped Allen key	1
5000222	Twist drill Ø 3.5 mm	1
5000224	Twist drill Ø 4.5 mm	1
5006010	Drill guide Ø 3.5 mm L 110 mm	1
5006.04	Ø 4.0 mm adapter t-wrench	1
5006.05	Ø 5.0 mm adapter t-wrench	1
5006.06	Ø 6.0 mm adapter t-wrench	1
5006564	ER rod template	1
5006560ER	ClickIt ER sterilization trax	1

Screw and wire range available by scanning the following QR code



Sterilization

The components of the ClickIt ER system are supplied sterile.

We indicate sterilizing the instruments in saturated steam autoclave according to the following parameters:

TEMPERATURE: 134°

EXPOSURE: 5 min

* Implantable and sterile devices are covered by the CE-1936 marking, while non-sterile reusable instruments are covered by the CE marking



HEADQUARTERS

Via P. Gobetti, 56r - 16145 Genova
Tel. +39 010 30801 - Fax +39 010 3080210
www.mikai.it - servizio.clienti@mikai.it

MANUFACTURING FACILITY

Via Canestrello, 2 - 36050 Monteviale (VI)
Tel. +39 0444 950100 - Fax +39 0444 950133
www.mikai.it - servizio.clienti@mikai.it