

FEP

MODULAR EXTERNAL FIXATOR

Lightweight
Modular
Sterile





Indications

Open, closed and intra-articular fractures.
Nonunions, limb lengthenings, deformity corrections and axial corrections.

FEATURES

Sterile kit

Available in three sizes


Versatile

CLINICAL ADVANTAGES


Ready for use

Wide range of application


Allows independent corrections on all planes



Kit FEP



Small FEP kit - 5000060MST		Medium FEP kit - 5000001MST		Large FEP kit - 5000030MST	
Description	Qty	Description	Qty	Description	Qty
Small fixator body	1	Medium fixator body	1	Large fixator body	1
Self-drilling cortical screw Ø 4.0 mm L 120 mm	4	Self-drilling cortical screw Ø 5.0 mm L 140 mm	6	Self-drilling cortical screw Ø 6.0 mm L 180 mm	6
Single clamp for FEP Small and Medium	4	Single clamp for FEP Small and Medium	6	Single clamp for FEP Large	6



FEP kit allows independent corrections in all 3 planes. The joint unit allows a procurvatum-recurvatum correction. The compression-distraction unit allows the correction of rotational defects. The micrometric regulators allow a varus-valgus correction. Any corrections can be made independently without losing the previous correction.

The elasticity of Aluminium alloy 7012 favors the reparative processes of the fracture site while ensuring an optimal stability.

Implantable and sterile devices are covered by the CE-1936 marking.



HEADQUARTERS

Via P. Gobetti, 56r - 16145 Genova
Tel. +39 010 30801 - Fax +39 010 3080210
www.mikai.it - servizio.clienti@mikai.it

MANUFACTURING FACILITY

Via Canestrello, 2 - 36050 Monteviale (VI)
Tel. +39 0444 950100 - Fax +39 0444 950133
www.mikai.it - servizio.clienti@mikai.it